


Product Datasheet

Human Fetuin A/AHSG ELISA Kit (orb50196)

Description	Human Fetuin A/AHSG ELISA Kit
Reactivity	Human
Range	0.78 ng/ml - 50 ng/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: A19-V367
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	39325 MW
Uniprot ID	P02765
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:100000, concentration in serum and plasma is 219-805 ug/ml.
Expiration Date	Please enquire.



Human
Fetuin A
PicoKine
ELISA Kit
standa...