


Product Datasheet

Mouse Galectin-1 ELISA Kit (orb50197)

Description	Mouse Galectin-1 ELISA Kit
Reactivity	Mouse
Range	156 pg/ml - 10,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: A2-E135
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	14866 MW
Uniprot ID	P16045
Sensitivity	5 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum and plasma is 1.6-2 ng/ml.
Expiration Date	Please enquire.



Mouse
Galectin-1
PicoKine
ELISA Kit
stan...