


Product Datasheet

Mouse IL-7/Interleukin-7 ELISA Kit (orb50198)

Description	Mouse IL-7/Interleukin-7 ELISA Kit
Reactivity	Mouse
Range	15.6 pg/ml - 1,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: E26-I154
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (EDTA).
MW	17727 MW
Uniprot ID	P10168
Sensitivity	1 pg/ml
Expiration Date	Please enquire.



Mouse IL-7
PicoKine
ELISA Kit
standard
C...