


Product Datasheet

Human Prostate-specific antigen ELISA kit (PSA/KLK3) (orb50199)

Description	Human Prostate-specific antigen ELISA kit (PSA/KLK3)
Reactivity	Human
Range	312 pg/ml - 20,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: A18-P261
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	28741 MW
Uniprot ID	P07288
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum and plasma is 0-4.78 ng/ml.
Expiration Date	Please enquire.



Human
Kallikrein 3
PicoKine
ELISA Kit
st...