


Product Datasheet

Human RBP4/Retinol Binding Protein 4 ELISA Kit (orb50200)

Description	Human RBP4/Retinol Binding Protein 4 ELISA Kit
Reactivity	Human
Range	312 pg/ml - 20,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: E19-L201
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum, plasma(heparin, EDTA) and urine.
MW	23010 MW
Uniprot ID	P02753
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:20000, concentration in serum and plasma is around 28 ug/ml.
Expiration Date	Please enquire.



Human
RBP4
PicoKine
ELISA Kit
standard
C...