


Product Datasheet

Mouse TrkB ELISA Kit (orb50201)

Description	Mouse TrkB ELISA Kit
Reactivity	Mouse
Range	93.8 pg/ml - 6,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: C32-H429
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA)
MW	92133 MW
Uniprot ID	P15209
Sensitivity	5 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum and plasma is 2-5 ng/ml.
Expiration Date	Please enquire.



Mouse TrkB
PicoKine
ELISA Kit
standard c...