


Product Datasheet

Rat TrkA ELISA Kit (orb50204)

Description	Rat TrkA ELISA Kit
Reactivity	Rat
Range	78 pg/ml - 5,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: A33-P418
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants and cell lysates.
MW	87868 MW
Uniprot ID	P35739
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum is around 1ng/ml
Expiration Date	Please enquire.



Rat TrkA
PicoKine
ELISA Kit
standard
cur...