


Product Datasheet

Human Laminin ELISA Kit (orb50206)

Description	Human Laminin ELISA Kit
Reactivity	Human
Range	156 pg/ml - 10,000 pg/ml
Tested Applications	ELISA
Immunogen	Target protein isolated from human fibroblasts
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates and serum.
MW	337084 MW
Uniprot ID	A4D0S4, P24043, Q16363, Q13751, Q13753, P25391, Q16787, P11047, Q9Y6N6, P07942, P55268, O15230
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:10, concentration in serum is 5-9 ng/ml.
Expiration Date	Please enquire.



Human Laminin ELISA Kit standard curve