


Product Datasheet

Mouse E-Cadherin/CDH1/Cadherin-1 ELISA Kit (orb50209)

Description	Mouse E-Cadherin/CDH1/Cadherin-1 ELISA Kit
Reactivity	Mouse
Range	156 pg/ml - 10,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: D157-V709
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin).
MW	98256 MW
Uniprot ID	P09803
Sensitivity	15 pg/ml
Dilution Range	Dilution ratio of 1:200, concentration in serum and plasma is 262-339 ng/ml.
Expiration Date	Please enquire.



Mouse E-Cadherin PicoKine ELISA Kit stan...