

Product Datasheet

**beta Catenin (phospho-Thr41/Ser45) antibody
(Biotin) (orb502092)**

Description

beta Catenin (phospho-Thr41/Ser45)

Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Biotin
Tested Applications	ELISA, ICC, IHC-Fr, IHC-P
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human Beta-Catenin around the phosphorylation site of Thr41/Ser45
Target	Phospho-Beta Catenin (Thr41 + Ser45)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	86kDa
Uniprot ID	P35222
Dilution Range	IHC-P=1:100-500,IHC-F=1:100-500,ICC=1:100,ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.