


Product Datasheet

Human CD56/NCAM-1 ELISA Kit (orb50211)

Description	Human CD56/NCAM-1 ELISA Kit
Reactivity	Human
Range	156 pg/ml - 10,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: L20-G718
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	94574 MW
Uniprot ID	P13591
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum and plasma is around 200 ng/ml.
Expiration Date	Please enquire.



Human
CD56/NCAM-
1 PicoKine
ELISA Kit
sta...