


Product Datasheet

Human CXCL10/IP-10 ELISA Kit (orb50212)

Description	Human CXCL10/IP-10 ELISA Kit
Reactivity	Human
Range	31.2 pg/ml - 2,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: V22-P98
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum, plasma (heparin, EDTA) , saliva and urine.
MW	10881 MW
Uniprot ID	P02778
Sensitivity	1 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum and plasma is 80-133 pg/ml.
Expiration Date	Please enquire.



Human
CXCL10/IP-
10 PicoKine
ELISA Kit
st...