


Product Datasheet

Human Mer ELISA Kit (orb50213)

Description	Human Mer ELISA Kit
Reactivity	Human
Range	156 pg/ml - 10,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: R26-A499
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates and serum.
MW	110249 MW
Uniprot ID	Q12866
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum is 896-3260 pg/ml.
Expiration Date	Please enquire.



Human Mer
PicoKine
ELISA Kit
standard
cu...