

Product Datasheet

Crk p38 (phospho Y221) antibody (Biotin) [Out of stock] (orb502139)

Description

Rabbit polyclonal antibody against

conjugated to Biotin

Species/Host

Rabbit

Conjugation

Biotin

Tested

ELISA, IF, IHC-Fr, IHC-P, WB

Applications
Immunogen

KLH conjugated Synthesised phosphopeptide derived from human Crk p38 around the phosphorylation site of Tyr221 [GP(p-Y)AQ]

Preservatives

0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.

Form/Appearance

Lyophilized or Liquid

Concentration

- 100 µg (in 100 µl): 1 mg/ml

Storage

Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

MW

38kDa

Uniprot ID
Q96HJ0, Q64010, P46108, Q63768, Q96GA9
Dilution Range

WB: 1:50-200, ELISA: 1:100-1000, IHC-P: 1:50-200, IF: 1:50-200

Expiration Date

12 months from date of receipt.