


Product Datasheet

Mouse LOX-1/OLR1 ELISA Kit (orb50214)

Description	Mouse LOX-1/OLR1 ELISA Kit
Reactivity	Mouse
Range	31.2 pg/ml - 2,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: R60-I363
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	41643 MW
Uniprot ID	Q9EQ09
Sensitivity	1 pg/ml
Expiration Date	Please enquire.



Mouse LOX-1/OLR1 PicoKine ELISA Kit stan...