

Product Datasheet

JAK2 (phospho-Tyr221) antibody (Biotin) (orb502163)

Description	JAK2 (phospho-Tyr221) antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human, Zebrafish
Conjugation	Biotin
Tested Applications	ELISA, ICC, IHC-Fr, IHC-P
Immunogen	KLH conjugated synthesised phosphopeptide derived from human JNK2 around the phosphorylation site of Tyr221
Target	phospho-JAK2 (Tyr221)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	O60674
Dilution Range	IHC-P=1:100-500,IHC-F=1:100-500,ICC=1:100,ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.