


Product Datasheet

Human PECAM-1/CD31 ELISA Kit (orb50218)

Description	Human PECAM-1/CD31 ELISA Kit
Reactivity	Human
Range	156 pg/ml - 10,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: CHO; Immunogen sequence: Q28-K601
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	82536 MW
Uniprot ID	P16284
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:10, concentration in serum and plasma is around 10 ng/ml.
Expiration Date	Please enquire.



Human
PECAM-
1/CD31
PicoKine
ELISA Kit
st...