

Product Datasheet

VEGFR2 (phospho-Tyr1059) antibody (Biotin) (orb502183)

Description	VEGFR2 (phospho-Tyr1059) antibody (Biotin)
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Biotin
Tested Applications	ELISA, IHC-P, WB
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human VEGF receptor 2 around the phosphorylation site of Tyr1059
Target	phospho-VEGF receptor 2 (Tyr1059)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P35968
Dilution Range	WB=1:500-2000,IHC-P=1:100-500,ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.