

Product Datasheet

JUP (phospho-Tyr550) antibody (Biotin) (orb502391)

Description	JUP (phospho-Tyr550) antibody (Biotin)
Species/Host	Rabbit
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human JUP around the phosphorylation site of Tyr550
Target	phospho-gamma Catenin (Tyr550)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P14923
Expiration Date	12 months from date of receipt.