

## Product Datasheet

**phospho-Tau proteinSer202 antibody (Biotin) [Out of stock] (orb502519)**

**Description**

Rabbit polyclonal antibody against

conjugated to Biotin

**Species/Host**

Rabbit

**Conjugation**

Biotin

**Tested**

ELISA, ICC, IF, IHC-Fr, IHC-P, WB

**Applications**
**Immunogen**

KLH conjugated Synthesised phosphopeptide derived from human Tau protein isoform 2 around the phosphorylation site of Ser202

**Preservatives**

0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.

**Form/Appearance**

Lyophilized or Liquid

**Concentration**

- 100 µg (in 100 µl): 1 mg/ml

**Storage**

Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.

**Note**

For research use only

**Isotype**

IgG

**Clonality**

Polyclonal

**MW**

52/79kDa

**Uniprot ID**
**P10637, P29172, P10636, P19332**
**Dilution Range**

WB: 1:50-200, ELISA: 1:100-1000, IHC-P: 1:50-200, ICC: 1:50-200, IF: 1:50-200

**Expiration Date**

12 months from date of receipt.