

Product Datasheet

Phospho-Dab1 Tyr198 antibody (Biotin) [Out of stock] (orb502537)

Description

Rabbit polyclonal antibody against

to Biotin

Species/Host

Rabbit

Reactivity

Human, Rat

Conjugation

Biotin

Tested

WB

Applications
Immunogen

KLH conjugated synthesised phosphopeptide derived from human Dab1 around the phosphorylation site of Tyr198 V(p-Y)QY

Preservatives

0.01M TBS (pH7.4) with 1% rAlbumin, 0.02% Proclin300 and 50% Glycerol.

Form/Appearance

Lyophilized or Liquid

Concentration

- 100 µg (in 100 µl): 1 mg/ml

Storage

Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

MW

65kDa

Uniprot ID
075553
Dilution Range

WB: 1:50-200, ELISA: 1:100-1000, IHC-P: 1:50-200, ICC: 1:50-200, IF: 1:50-200

Expiration Date

12 months from date of receipt.