

Product Datasheet

RAC1 antibody (Biotin) (orb502548)

Description	RAC1 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA, ICC
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human RAC1 around the phosphorylation site of Ser71
Target	phospho-RAC1 + RAC2 (Ser71)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P63000
Dilution Range	ICC=1:100,ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.