

Product Datasheet

phospho-GluA4 Ser862 antibody (FITC) [Out of stock] (orb502607)

Description	Rabbit polyclonal antibody to Anti-phospho-GLUR4 (Ser862)/FITC
Species/Host	Rabbit
Conjugation	FITC
Immunogen	KLH conjugated synthesised phosphopeptide derived from human GluA4 around the phosphorylation site of Ser862
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Lyophilized or Liquid
Concentration	1mg/ml
Storage	Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	98kDa
Uniprot ID	P19493, P48058, Q9Z2W8
Expiration Date	12 months from date of receipt.