

Product Datasheet

Beta-Catenin (phospho-Ser552) antibody (FITC) (orb502722)

Description	Beta-Catenin (phospho-Ser552) antibody (FITC)
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	FITC
Tested Applications	FC, ICC, IF
Immunogen	KLH conjugated synthesised phosphopeptide derived from human beta Catenin around the phosphorylation site of Ser552
Target	Phospho-Beta-Catenin (Ser552)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P35222
Dilution Range	Flow-Cyt=1µg/Test,IF=1:100-500,ICC=1:100-500
Expiration Date	12 months from date of receipt.