

Product Datasheet

SMC1 (phospho-Ser966) antibody (FITC) [Out of stock] (orb502750)

Description

Rabbit polyclonal antibody against

to FITC

Species/Host

Rabbit

Conjugation

FITC

Tested

ICC, IF

Applications
Immunogen

KLH conjugated synthesised phosphopeptide derived from human SMC1 around the phosphorylation site of Ser966.

Preservatives

0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.

Form/Appearance

Lyophilized or Liquid

Concentration

- 100 µg (in 100 µl): 1 mg/ml

Storage

Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

MW

136kDa

Uniprot ID
Q9CU62, Q8NDV3, Q14683, Q920F6
Dilution Range

ICC: 1:50-200, IF: 1:50-200

Expiration Date

12 months from date of receipt.