

Product Datasheet

Phospho-c-Fos Ser374 antibody (FITC) (orb502838)

Description

Rabbit polyclonal antibody against

FITC

Species/Host

Rabbit

Reactivity

Mouse, Rat

Conjugation

FITC

Tested

FC

Applications
Immunogen

KLH conjugated Synthesised phosphopeptide derived from human c-Fos around the phosphorylation site of Ser374

Target

FOS

Preservatives

0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.

Form/Appearance

Lyophilized or Liquid

Concentration

- 100 µg (in 100 µl): 1 mg/ml

Storage

Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

MW

41kDa

Uniprot ID
P01100
Dilution Range

Flow-Cyt=,IF=1:100-500,ICC=1:100-500

Expiration Date

12 months from date of receipt.