

Product Datasheet

**phospho-Tau proteinThr522 antibody (FITC)
(orb502921)**

Description

Rabbit polyclonal antibody against

conjugated to FITC

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

FITC

Tested

FC

Applications
Immunogen

KLH conjugated synthesised phosphopeptide derived from human Tau protein around the phosphorylation site of Thr522

Target

MAPT

Preservatives

0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.

Form/Appearance

Lyophilized or Liquid

Concentration

- 100 µg (in 100 µl): 1 mg/ml

Storage

Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

MW

52/79kDa

Uniprot ID
P10636
Dilution Range

Flow-Cyt=,IF=1:100-500,ICC=1:100-500

Expiration Date

12 months from date of receipt.