

Product Datasheet

Phospho-e2 alphaSer51 antibody (FITC) (orb502926)

Description

Rabbit polyclonal antibody against

FITC

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

FITC

Tested

IF

Applications
Immunogen

KLH conjugated synthesised phosphopeptide derived from human eIF2 alpha around the phosphorylation site of Ser51 [EL(p-S)RR]

Target

EIF2S1

Preservatives

0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.

Form/Appearance

Lyophilized or Liquid

Concentration

- 100 µg (in 100 µl): 1 mg/ml

Storage

Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

MW

36kDa

Uniprot ID
P05198
Expiration Date

12 months from date of receipt.