

Product Datasheet

PI3K (phospho-Tyr524) antibody (HRP) (orb503041)

Description

PI3K (phospho-Tyr524) antibody (HRP)

Species/Host

Rabbit

Reactivity

Mouse

Conjugation

HRP

Tested

ELISA, IHC-Fr, IHC-P, WB

Applications
Immunogen

KLH conjugated Synthesised phosphopeptide derived from human PIK3CD around the phosphorylation site of Tyr524

Target

phospho-PI 3 Kinase p110 delta (Tyr524)

Preservatives

0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.

Form/Appearance

Liquid

Concentration

1mg/ml

Storage

Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

MW

115kDa

Uniprot ID
[O00329](https://www.uniprot.org/entry/O00329)
Dilution Range

WB=1:500-2000,IHC-P=1:100-500,IHC-F=1:100-500,ELISA=1:5000-10000

Expiration Date

12 months from date of receipt.