

Product Datasheet

MAPK8 (phospho-Thr183+Tyr185) antibody (HRP) (orb503086)

Description	MAPK8 (phospho-Thr183+Tyr185) antibody (HRP)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	HRP
Tested Applications	ELISA, IHC-Fr, IHC-P
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human MAPK9 around the phosphorylation site of Thr183+Tyr185
Target	phospho-MAPK8 (Thr183 + Tyr185)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P45983
Dilution Range	IHC-P=1:100-500,IHC-F=1:100- 500,ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.