

Product Datasheet

VEGF receptor 2 (phospho-Tyr996) antibody (HRP) (orb503266)

Description

Rabbit polyclonal antibody to Anti-

(Tyr996)/HRP

Species/Host

Rabbit

Conjugation

HRP

Tested

WB

Applications
Immunogen

KLH conjugated Synthesised phosphopeptide derived from human VEGF receptor 2 around the phosphorylation site of Tyr996

Target

KDR

Preservatives

0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.

Form/Appearance

Lyophilized or Liquid

Concentration

1mg/ml

Storage

Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

MW

147kDa

Uniprot ID
P35968
Expiration Date

12 months from date of receipt.