

Product Datasheet

**PLC gamma 2 (phospho-Tyr1217) antibody (HRP)
(orb503333)**

Description

Rabbit polyclonal antibody against

conjugated to HRP

Species/Host

Rabbit

Reactivity

Human

Conjugation

HRP

Tested

ELISA, IHC-Fr, IHC-P, WB

Applications
Immunogen

KLH conjugated Synthesised phosphopeptide derived from human PLC gamma 2 around the phosphorylation site of Tyr1217

Target

PLCG2

Preservatives

0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.

Form/Appearance

Lyophilized or Liquid

Concentration

- 100 µg (in 100 µl): 1 mg/ml

Storage

Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

MW

139kDa

Uniprot ID
P16885
Dilution Range

WB=1:500-2000,IHC-P=1:100-500,IHC-F=1:100-500,ELISA=1:5000-10000

Expiration Date

12 months from date of receipt.