

Product Datasheet

HSP27 (phospho-Ser82) antibody (HRP) (orb503531)

Description	HSP27 (phospho-Ser82) antibody (HRP)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	HRP
Tested Applications	ELISA, IHC-P
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human HSP27 around the phosphorylation site of Ser82
Target	Phospho-HSP27 (Ser82)
Preservatives	0.01M PBS, 0.5%rAlbumin, 0.03% Proclin300
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P04792
Dilution Range	IHC-P=1:100-500,ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.