

## Product Datasheet

### EPB41L2 Immunizing peptide (orb50354)

|                          |                                                                              |
|--------------------------|------------------------------------------------------------------------------|
| <b>Description</b>       | EPB41L2 Immunizing peptide                                                   |
| <b>Preservatives</b>     | 100µg of Lyophilised peptide, formulation<br>De-ionized (DI) water.          |
| <b>Storage</b>           | Shipped at ambient temperature, store at -<br>20°C                           |
| <b>Note</b>              | For research use only                                                        |
| <b>Application notes</b> | This peptide was used as immunogen to<br>produce orb50354 (EPB41L2 antibody) |
| <b>Protein Sequence</b>  | DYDPEEHGSID                                                                  |
| <b>Expiration Date</b>   | 12 months from date of receipt.                                              |