

Product Datasheet

phospho-alpha 1 Sodium Potassium ATPase Tyr260 antibody (HRP) (orb503577)

Description	phospho-alpha 1 Sodium Potassium ATPase Tyr260 antibody (HRP)
Species/Host	Rabbit
Reactivity	Human
Conjugation	HRP
Tested Applications	ELISA, ICC, IHC-Fr, WB
Immunogen	KLH conjugated synthesised phosphopeptide derived from human alpha 1 Sodium Potassium ATPase around the phosphorylation site of Tyr260
Target	phospho-ATP1A1 (Tyr260)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P05023
Dilution Range	WB=1:500-2000,IHC-F=1:100-500,ICC=1:100,ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.