

Product Datasheet

**14-3-3 beta+zeta (phospho-Ser186) antibody
(HRP) (orb503619)**

Description

14-3-3 beta+zeta (phospho-Ser186)

Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	HRP
Tested Applications	ELISA, IHC-Fr, IHC-P
Immunogen	KLH conjugated synthesised phosphopeptide derived from human 14-3-3 beta + zeta around the phosphorylation site of Ser186 / Ser184
Target	phospho-14-3-3 beta + zeta (Ser186 / Ser184)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P31946
Dilution Range	IHC-P=1:100-500,IHC-F=1:100-500,ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.