

Product Datasheet

GRIN2B (phospho-Tyr1472) antibody (HRP) [Out of stock] (orb503693)

Description

Rabbit polyclonal antibody against

conjugated to HRP

Species/Host

Rabbit

Conjugation

HRP

Tested

ELISA, ICC, IHC-Fr, IHC-P, WB

Applications
Immunogen

KLH conjugated Synthesised phosphopeptide derived from Rat NMDAR2B around the phosphorylation site of Tyr1472

Preservatives

0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.

Form/Appearance

Lyophilized or Liquid

Concentration

- 100 µg (in 100 µl): 1 mg/ml

Storage

Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

MW

180kDa

Uniprot ID
Q00960, Q01097, Q13224
Dilution Range

WB: 1:50-200, ELISA: 1:100-1000, IHC-P: 1:50-200, ICC: 1:50-200

Expiration Date

12 months from date of receipt.