

Product Datasheet

C-Myc (phospho-Ser62) antibody (HRP) (orb503728)

Description	C-Myc (phospho-Ser62) antibody (HRP)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	HRP
Tested Applications	ELISA, IHC-P, WB
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human C-Myc around the phosphorylation site of Ser62
Target	phospho-C-Myc (Ser62)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P01106
Dilution Range	WB=1:500-2000,IHC-P=1:100-500,ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.