

Product Datasheet

IRS1 (phospho-Tyr1179) antibody (HRP) (orb503742)

Description	IRS1 (phospho-Tyr1179) antibody (HRP)
Species/Host	Rabbit
Reactivity	Mouse, Zebrafish
Conjugation	HRP
Tested Applications	ELISA, IHC-Fr, IHC-P
Immunogen	KLH conjugated synthesised phosphopeptide derived from human IRS1 around the phosphorylation site of Ser1179
Target	phospho-IRS1 (Tyr1179)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P35568
Dilution Range	IHC-P=1:100-500,IHC-F=1:100-500,ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.