

Product Datasheet

ZIP Kinase (phospho-Thr265) antibody (HRP) (orb503744)

Description

ZIP Kinase (phospho-Thr265) antibody

Species/Host	Rabbit
Reactivity	Mouse, Rat
Conjugation	HRP
Tested Applications	ELISA, ICC, IHC-Fr, IHC-P, WB
Immunogen	KLH conjugated synthesised phosphopeptide derived from human ZIP Kinase around the phosphorylation site of Thr265
Target	phospho-ZIP Kinase (Thr265)
Preservatives	Preservative: 15mM Sodium Azide, Constituents: 1% rAlbumin, 0.01M PBS, pH 7.4
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	O43293
Dilution Range	WB=1:500-2000,IHC-P=1:100-500,IHC-F=1:100-500,ICC=1:100-500,ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.