

Product Datasheet

phospho-FANCD2Ser1404 antibody (HRP) [Out of stock] (orb503894)

Description

Rabbit polyclonal antibody against

conjugated to HRP

Species/Host

Rabbit

Conjugation

HRP

Tested

ELISA, ICC, IHC-Fr, IHC-P

Applications
Immunogen

KLH conjugated Synthesised phosphopeptide derived from human FANCD2 around the phosphorylation site of Ser1404

Preservatives

0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.

Form/Appearance

Lyophilized or Liquid

Concentration

- 100 µg (in 100 µl): 1 mg/ml

Storage

Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

MW

162kDa

Uniprot ID
Q80V62, Q6IV68, Q9BXW9
Dilution Range

ELISA: 1:100-1000, IHC-P: 1:50-200, ICC: 1:50-200

Expiration Date

12 months from date of receipt.