

Product Datasheet

phospho-Ephexin-1 Tyr179 antibody (HRP) (orb503906)

Description	phospho-Ephexin-1 Tyr179 antibody (HRP)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	HRP
Tested Applications	ELISA, WB
Immunogen	KLH conjugated synthesised phosphopeptide derived from human Ephexin-1 around the phosphorylation site of Tyr179
Target	phospho-Ephexin-1 (Tyr179)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	82kDa
Uniprot ID	Q8N5V2
Dilution Range	WB=1:500-2000,ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.