

Product Datasheet

Phospho-EEF2 Thr56 + Thr58 antibody (HRP) (orb503909)

Description	Phospho-EEF2 Thr56 + Thr58 antibody (HRP)
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	HRP
Tested Applications	ELISA, WB
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human EEF2 around the phosphorylation site of Thr56+Thr58
Target	Phospho-EEF2 (Thr56 + Thr58)
Preservatives	Preservative: 15mM Sodium Azide, Constituents: 1% rAlbumin, 0.01M PBS, pH 7.4
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P13639
Dilution Range	WB=1:500-2000,ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.