

Product Datasheet

PDK1 (phospho-Tyr373/376) antibody (HRP) (orb504028)

Description	PDK1 (phospho-Tyr373/376) antibody (HRP)
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	HRP
Tested Applications	ELISA, WB
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human PDPK1 around the phosphorylation site of (Y373+Y376)
Target	Phospho-PDPK1 (Tyr373 + Tyr376)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	O15530
Dilution Range	WB=1:500-2000,ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.