

Product Datasheet

**p53 (phospho-Ser6) antibody (HRP) [Out of stock]
(orb504198)**

Description

Rabbit polyclonal antibody against

HRP

Species/Host

Rabbit

Conjugation

HRP

Tested

ELISA

Applications
Immunogen

KLH conjugated Synthesised phosphopeptide derived from human P53 around the phosphorylation site of Ser6

Preservatives

0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.

Form/Appearence

Lyophilized or Liquid

Concentration

- 100 µg (in 100 µl): 1 mg/ml

Storage

Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

MW

43kDa

Uniprot ID

P02340, **P04637**, **P10361**, **P67939**

Dilution Range

ELISA: 1:100-1000

Expiration Date

12 months from date of receipt.