

Product Datasheet

ELk1 (phospho-Ser389) antibody (HRP) (orb504285)

Description	ELk1 (phospho-Ser389) antibody (HRP)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	HRP
Tested Applications	ELISA, IHC-Fr, WB
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human ELk1 around the phosphorylation site of Ser389
Target	Phospho-ELk1 (Ser389)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P19419
Dilution Range	WB=1:500-2000,IHC-F=1:100-500,ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.