

## Product Datasheet

**PIK3C3 (phospho-Ser164) antibody (HRP) [Out of stock] (orb504375)**

**Description**

Rabbit polyclonal antibody against

conjugated to HRP

**Species/Host**

Rabbit

**Conjugation**

HRP

**Tested**

ELISA, IHC-Fr, IHC-P, WB

**Applications**
**Immunogen**

KLH conjugated Synthesised phosphopeptide derived from human PIK3C3 around the phosphorylation site of Ser164

**Preservatives**

0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.

**Form/Appearance**

Lyophilized or Liquid

**Concentration**

- 100 µg (in 100 µl): 1 mg/ml

**Storage**

Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.

**Note**

For research use only

**Isotype**

IgG

**Clonality**

Polyclonal

**MW**

98kDa

**Uniprot ID**
**O88763, Q8NEB9, Q6PF93**
**Dilution Range**

WB: 1:50-200, ELISA: 1:100-1000, IHC-P: 1:50-200

**Expiration Date**

12 months from date of receipt.