

Product Datasheet

Phospho-Dab1 Tyr198 antibody (HRP) [Out of stock] (orb504475)

Description

Rabbit polyclonal antibody against

	HRP
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	HRP
Tested Applications	WB
Immunogen	KLH conjugated synthesised phosphopeptide derived from human Dab1 around the phosphorylation site of Tyr198 V(p-Y)QY
Preservatives	0.01M TBS (pH7.4) with 1% rAlbumin, 0.02% Proclin300 and 50% Glycerol.
Form/Appearance	Lyophilized or Liquid
Concentration	- 100 µg (in 100 µl): 1 mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	65kDa
Uniprot ID	075553
Dilution Range	WB: 1:50-200, ELISA: 1:100-1000, IHC-P: 1:50-200, ICC: 1:50-200
Expiration Date	12 months from date of receipt.