

Product Datasheet

CDK5 (phospho-Tyr15) antibody (PE) (orb504539)

Description	Rabbit polyclonal antibody against CDK5 (phospho-Tyr15) conjugated to PE
Species/Host	Rabbit
Reactivity	Mouse, Rat
Conjugation	PE
Tested Applications	IF
Immunogen	KLH conjugated synthesised phosphopeptide derived from human CDK5 around the phosphorylation site of Tyr15
Target	CDK5
Preservatives	0.01M TBS(pH7.4) with 0.03% Proclin300.
Form/Appearance	Lyophilized or Liquid
Concentration	- 100 µg (in 100 µl): 1 mg/ml
Storage	Store at 2-8 °C for one year. Protect from light.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	32kDa
Uniprot ID	Q00535
Expiration Date	12 months from date of receipt.