

Product Datasheet

APBB1 (phospho-Ser175) antibody (PE) [Out of stock] (orb504768)

Description

Rabbit polyclonal antibody to Anti-

Species/Host	Rabbit
Conjugation	PE
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human APBB1 around the phosphorylation site of Ser175 [LS(p-S)PP]
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Lyophilized or Liquid
Concentration	1mg/ml
Storage	Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	77kDa
Uniprot ID	P46933, Q9QXJ1, O00213
Dilution Range	IF=1:50-200 not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user.
Expiration Date	12 months from date of receipt.